



## BlasTaq™ 2X PCR MasterMix

### Cat. No. G895

Store at -20°C.

### Product Description

**BlasTaq™ 2X PCR MasterMix** is a ready-to-use MasterMix containing **abm's** BlasTaq™ DNA Polymerase in a uniquely-formulated buffer with gel loading dye. This strategically-engineered, next generation Taq Polymerase provides **rapid extension rates and robust performance**. With specialized reaction conditions, this polymerase provides increased processivity, yields, and sensitivity, while shortening reaction times by up to 70%, compared to wild-type Taq DNA polymerase. BlasTaq™ has 5'-3' polymerase and 5'-3' exonuclease activities, lacks 3'-5' exonuclease activity, and produces 3'-dA-tailed amplicons. PCR products made with BlasTaq™ can be used with TA cloning vectors.

Product Component	Quantity	Cat. No.
BlasTaq™ 2X PCR MasterMix <sup>1</sup>	800 rxn (10.0 ml)	G895

<sup>1</sup> Buffer contains 1.5 mM Mg<sup>2+</sup>.

### Protocol

1. Mix individual components before use and assemble reaction on ice.

Component	Volume
2X BlasTaq™ PCR MasterMix	12.5 µl
Forward Primer (10 µM)	1 µl
Reverse Primer (10 µM)	1 µl
Template DNA	Variable (100 ng genomic DNA)
Nuclease-free H <sub>2</sub> O	Up to 25 µl

2. Gently mix the reaction components and briefly centrifuge. Run thermocycling conditions for standard PCR:

Step	Temperature	Duration
Initial Denaturation <sup>2</sup>	95°C	3 min
25 – 35 Cycles	95°C	15 sec
	60°C <sup>3</sup>	15 sec
	72°C	15 sec/kb
Final Extension	72°C	1 min

<sup>2</sup> For most applications, an initial 3 minute denaturation step at 95°C is sufficient. Increase to 5 minutes for high-GC or difficult templates.

<sup>3</sup> BlasTaq™'s PCR buffer allows for primer annealing at 60°C for most primers and adjust only if needed.

3. After PCR, maintain the reaction at 4°C or store at -20°C until use.
4. Analyze the amplification products by agarose gel electrophoresis.
5. Visualize by ethidium bromide or SafeView™ (Cat No. **G108**) staining.

### General Notes

- Specialized buffer for higher yields, sensitivity, and specificity compared to wild-type Taq polymerase.
- Decrease reaction times by 70% using specialized protocol.
- For optimal efficiency, use a 25 µl reaction volume.