

Purpose:

To prepare **LAD2/LADR cells** for the **Mast Granulation Assay (G7800)** by sensitization with human IgE (Y70000) and crosslinking with anti-IgE (Y79000).

Materials:

Cat. No.	Product
T8157/ T8156	LAD2/ LADR cells
Y70000	Human IgE
Y79000	Human anti-IgE

Additionally you will require:

- hSCF-free culture medium (*complete medium*)
- hSCF-free medium + 0.5% BSA (*washing solution*)
- 15 mL conical tubes
- 6-well plates
- A centrifuge

Protocol:

1. Cell preparation:

- Thaw and recover LAD2/LADR cells following standard culture guidelines.
- Count cells once they have recovered; use $1-10 \times 10^6$ cells per condition; the example below uses 2×10^6 cells/ condition.
- Removal of hSCF:
 - Collect $> 3 \times 10^6$ cells per condition in 15 mL tubes.
 - Centrifuge at 1500 rpm for 5 min and discard supernatant.
 - Resuspend in 2 mL and wash pellet 3 times with hSCF-free medium.
 - After final wash, resuspend in 1 mL hSCF-free medium and count cells.

2. Seeding:

- Seed 1×10^6 cells per well in 1 mL hSCF-free medium in 6-well plate:
 - Well 1: **“Activated” (IgE treatment)**
 - Well 2: **“Unactivated” (control)**
- Incubate 4 hours at 37 °C, 5% CO₂.

3. IgE sensitization:

- Add the IgE to the “Activated” well to reach a **final concentration of 2 µg/mL**. Leave the unactivated well in hSCF-free medium as a control.

- b. Incubate overnight at 37°C, 5% CO₂.

Note: If needed, pre-dilute IgE (Y70000) in hSCF-free medium before addition.

4. IgE crosslinking:

- a. Prepare washing solution: hSCF-free medium + 0.5% BSA.
- b. Wash activated cells 3 times with washing solution:
 - i. Centrifuge at 1500 rpm for 5 min each time.
 - ii. Resuspend pellet in 2 mL wash solution each wash.
- c. Resuspend the washed:
 - i. Activated cells in **1 mL hSCF-free medium** containing **anti-human IgE (Y79000) at a final concentration of 1 µg/mL** in 6-well plate.
- d. Incubate 1 hour at 37°C, 5% CO₂.

Note: If needed, pre-dilute anti-IgE (Y79000) in hSCF-free medium before addition.

Please be advised that the washing step here is very critical. Complete wash away of the un-bound human IgE is important to get the correct testing result.

5. Proceed to Mast Granulation Assay (G7800)

- a. After the 1-hour IgE crosslinking step, activated and control cells are ready for downstream degranulation assay.

Tips:

- Ensure complete removal of hSCF before IgE sensitization to avoid baseline activation.
- Include unactivated controls to monitor background degranulation