

Plasmid Amplification

Plasmid Amplification Protocol

abm supplies purified plasmid DNA in TE buffer (10 mM Tris, 1 mM EDTA, pH 8.0), optimized for efficient amplification through direct transformation into standard *E. coli* cloning strains such as DH5a or DH10B derivatives.

Transformation Protocol

The following protocol serves as a general guideline and may require optimization.

1. Thaw a vial of competent cells on ice for 10 min. Gently flick tube to thoroughly mix any settled cells. Aliquot 60 μ l of cells per transformation into sterile 1.5 ml tubes.
Note: once thawed, cells should not be re-frozen as transformation efficiency will be compromised.
2. Add 1 μ l of plasmid DNA to the cells and mix gently by flicking the tube 4-5 times. Do not vortex. Incubate reaction at 4°C for 30 min.
3. Heat shock the cells by placing the tube into a 42°C water bath for 45 s. Remove tubes immediately and place on ice for 1-2 min.
4. Add 150 μ l of room temperature SOC media to the mixture. Incubate tubes at 37°C with vigorous shaking (250 rpm) for 1 h.
5. Prepare several dilutions of the mixture and spread 50-100 μ l of each dilution onto appropriate antibiotic selection plates. Incubate at 37°C overnight.
6. Inoculate a single colony into 5-10 ml of LB broth with antibiotic selection. Incubate overnight at 200 rpm in a 37°C shaking incubator for 16-18 h.
7. Isolate plasmid DNA by standard miniprep protocol.

E. coli Antibiotic Selection Table

	Kanamycin	Ampicillin	Spectinomycin	Tetracycline	Chloramphenicol
Working Concentration	50 μ g/ml	100 μ g/ml	50 μ g/ml	12.5 μ g/ml	25 μ g/ml