


## LADR and LAD2 Functional Validation

**Product:** LADR Mast Cells ([Cat. No. T8156](#)) & LAD2 Mast Cells ([Cat. No. T8157](#))  
**Assay:** Mast Cell Granulation Assay ([Cat. No. G7800](#))

To evaluate the functional activity of mast cells, LADR and LAD2 cells were subjected to the Mast Granulation Assay following IgE sensitization with human IgE ([Cat. No. Y70000](#)) and IgE crosslinking ([Cat. No. Y79000](#)). Degranulation was quantified by measuring p-nitroaniline (pNA)\* release in the supernatant (released granules) and lysate (retained granules).

\*Note: p-Nitroaniline (pNA) is not histamine and is not a naturally occurring marker of histamine release. Instead, pNA is a chromogenic substrate used in laboratory assays to quantify tryptase activity, an enzyme co-released with histamine during mast cell degranulation. The amount of pNA generated is therefore an indirect but reliable indicator of mast cell activation.

### Results Summary



**Figure 1: IgE-mediated degranulation of LADR and LAD2 mast cells.** Only activated cells trigger a significant rise in pNA release, demonstrating that both cell lines retain robust, stimulus-dependent mast cell functionality.

## Data Analysis

Upon stimulation, both LADR and LAD2 cells demonstrated stimulus-dependent degranulation, with LADR showing a stronger overall response. Both cell lines maintain low baseline degranulation in the absence of stimulation.

LADR Observations	
Stimulated	<ul style="list-style-type: none"><li>• High pNA levels in the lysate, indicating a substantial store of granules.</li><li>• Significant increase in supernatant pNA, confirming robust granule release following activation.</li></ul>
Unstimulated	<ul style="list-style-type: none"><li>• Elevated lysate pNA levels consistent with stored granules.</li><li>• Minimal pNA release into the supernatant, confirming tight regulation of granulation without a trigger.</li></ul>

LAD2 Observations	
Stimulated	<ul style="list-style-type: none"><li>• Moderate lysate pNA levels with measurable release into the supernatant.</li><li>• Activation response is present, though lower in magnitude compared to LADR.</li></ul>
Unstimulated	<ul style="list-style-type: none"><li>• Lower lysate pNA levels compared to LADR, consistent with known cell-line differences.</li><li>• Supernatant pNA nearly absent, indicating minimal spontaneous release.</li></ul>

## Conclusion

Both LADR and LAD2 maintain stable, quiescent baseline behavior, with little to no spontaneous degranulation. Both cell lines respond to IgE crosslinking with measurable degranulation, with LADR exhibiting a stronger activation profile.