

qPCR Lentivirus Titer Kit

Cat. No. LV900

Store at -20°C.

Product Description

abm's qPCR Lentivirus Titer Kit is a one-step assay by qPCR which **does not require additional lysis or reverse transcription (RT) steps**, and can be completed under one hour. Designed to deliver **high sensitivity and specificity**, the kit ensures minimal non-specific background and better overall performance compared to similar kits on the market. ROX reference dye is provided separate from the MasterMix, making this kit universally compatible with most qPCR instruments.

| Product Component | Quantity | Part No. |
|----------------------------------|------------------|----------|
| BlasTaq™ 2X qPCR Titer MasterMix | 1.25 ml | P889-1 |
| Primer Mix | 100 rxn (200 µl) | LV900-A |
| Standard Control DNA | 50 µl | LV900-B |
| ROX Reference Dye | 15 µl | P101 |
| Nuclease-Free H ₂ O | 2 x 1.0 ml | P100 |

Protocol

The recommended amount of ROX Reference Dye to be added into the MasterMix may vary depending on the qPCR machine type:

- No ROX equipment: Not needed.
- Low ROX equipment: 1 µl/1.25 ml MasterMix.
- High ROX equipment: 11 µl/1.25 ml MasterMix.

1. **Sample Preparation:** For purified high titer viral samples, dilute the virus to 10⁷ IU/ml range with 1X PBS or DMEM. For low viral titer samples, collect viral supernatant for direct qPCR set up.
2. **Standard Control DNA Dilutions:** Perform five (5) 10-fold serial dilutions of the Standard Control DNA by diluting 5 µl DNA into 45 µl Nuclease-free H₂O in each step. Dilutions 1/100 to 1/100,000 will be used for generating the standard curve.

3. **Set-up:** All reactions are recommended to be set-up on ice in duplicates.

| Component | Volume |
|--------------------------------|--------|
| 2X qPCR MM | 10 µl |
| Primer Mix | 2 µl |
| Sample, NTC, or Standard DNAs | 2 µl |
| Nuclease-free H ₂ O | 6 µl |

4. qPCR cycling conditions:

| Step | Temperature | Duration | Cycle(s) |
|---------------------|-------------|----------|----------|
| Enzyme Activation | 95°C | 3 min | 1 |
| Denaturation | 95°C | 15 sec | 35 |
| Annealing/Extension | 60°C | 1 min | |

Data Analysis

Plot Ct value (**Y-axis**, linear scale) vs. Virus titer (**X-axis**, logarithmic scale). Generate a logarithmic regression using the four (4) Standard Control DNA dilutions to determine the unknown virus sample titer using $y = m \ln(x) + b$ from the trendline equation. The R² value should be > 0.95 to justify the proper assay setup. Note to include the dilution factor in the final calculation (i.e. if you diluted your purified viral sample 1/100 in Step 1, then the titer of the unknown sample should be multiplied by a factor of 100).

$$\text{Virus titer (IU/ml)} = e^{(Ct - b)/m}, \text{ where } m \text{ is the slope of the line and } b \text{ is the y-intercept.}$$

Example: trendline equation is $y = -1.349 \ln(x) + 40.898$; Ct of unknown sample = 16.98

$$\text{Virus titer (IU/ml)} = e^{(16.98 - 40.898)/-1.349} = 5.01 \times 10^7 \text{ IU/ml}$$

| Dilution | Virus Titer (IU/ml) |
|-----------|---------------------|
| 1/100 | 5 x 10 ⁸ |
| 1/1,000 | 5 x 10 ⁷ |
| 1/10,000 | 5 x 10 ⁶ |
| 1/100,000 | 5 x 10 ⁵ |

Download **abm's** calculation file from the product page on our website.