



2nd Generation Packaging Mix

Cat. No. LV003

Store at -20°C.

Product Description

abm's Second Generation Packaging Mix offers a high-safety, high-efficiency solution for lentivirus production. This system requires three plasmids for viral particle production: 1) Transfer Plasmid containing your gene of interest, 2) Packaging Plasmid containing structural proteins Gag, Pol, Tat, Rev, and 3) Envelope Plasmid containing transmembrane glycoprotein VSV-G. The ready-to-use 2nd Generation Packaging Mix combines optimized Packaging and Envelope Plasmids for high viral titers, enhanced biosafety and reliable gene delivery performance.

Packaging Mix Component	Quantity
pLenti-P2A	200 µl, 500 ng/µl
pLenti-P2B	

Additional Materials Required (not included)

Material	Recommended Product	Cat. No.
HEK293T Cell Line	293T Cells	LV010
Serum-Free DMEM	PriGrow III Medium	TM003
Transfection Reagent	DNAfectin™ Plus	G2500
Transduction Enhancer	ViralEntry™ Transduction Enhancer	G515

Protocol

Table 1: Reagent Requirements for Lentivirus Production

Vessel	Seeding Density	DNA Mix			Transfection Mix	
		Transfer Plasmid	Packaging Mix	Serum-Free DMEM	Transfection Reagent	Serum-Free DMEM
15 cm	12.5 x 10 ⁶	30 µg	60 µl	2.5 ml	160 µl	2.5 ml
10 cm	5.0 x 10 ⁶	10 µg	20 µl	1 ml	80 µl	1 ml
6-well	1.0 x 10 ⁶	2.5 µg	5 µl	100 µl	16 µl	100 µl

*Note: the specified amounts and protocol apply to the recommended products and may differ when using alternative products.

Day 1

- Seed an appropriate number of HEK293T cells into a desired vessel according to Table 1 using complete growth media, and incubate at 37°C with 5% CO₂ overnight.

Day 2

- Verify that cells have reached 70-80% confluence before proceeding with transfection.
- Based on your selected vessel, use Table 1 to prepare two solutions, the **DNA Mix** and **Transfection Mix** in 1.5 ml tubes. The DNA Mix includes: Transfer Plasmid, Packaging Mix and Serum-Free DMEM. The Transfection Mix includes: Transfection Reagent and Serum-Free DMEM. Incubate at room temperature for 5 min.
- Prepare the **Transfection Complex** by combining the DNA Mix and Transfection Mix together, and then incubate at room temperature for 20 min.
- Add Serum-Free DMEM to the Transfection Complex (800 µl for 6-well, 4 ml for 10 cm dish, 5 ml for 15 cm dish).
- Aspirate media from the vessel and gently add the Transfection Complex to the cells. Incubate at 37°C with 5% CO₂ for 5-8 h.
- Add complete growth media to the vessel (1 ml for 6-well, 4 ml for 10 cm dish, 5 ml for 15 cm dish) and incubate at 37°C with 5% CO₂ overnight.

Day 3

- Aspirate media from the vessel and add an appropriate amount of complete growth media. Incubate at 37°C with 5% CO₂ for 24 h.

Day 4

- Collect the supernatant from the vessel into a centrifuge tube and centrifuge at 1500 x g for 15 min at 4°C. Transfer the clarified supernatant to a fresh tube.
- Apply the clarified supernatant over a PES 0.45 µm sterile filter. Use lentivirus immediately or store at 4°C (short term) or -80°C (long term).
- Optional: A second harvest can be performed by adding an appropriate amount of complete growth media to the remaining cells from Day 4 Step 1 and incubating at 37°C with 5% CO₂ for an additional 24 h. The following day, perform Day 4 Step 1-2 procedure and combine the filtered supernatant with the first harvest.
- Recommended: Add **ViralEntry™ Transduction Enhancer** to the culture media when using the lentivirus to enhance downstream transduction efficiency.