



Certificate of Analysis

Product Description

Product Name	Immortalized Human Pancreatic Cells - Myc
Cat Number	T0158
Lot Number	VS7951
Source Organism	<i>Homo sapiens</i>
Source Organ	Pancreas
Viability	>80%
Freezing Medium	Cryopreservation Medium (TM024)
Date of Cryopreservation	June 24, 2014
QC Evaluation Medium	Prigrow II (TM002)

Test for Microbial Contaminants

	Test Method	Results
Bacteria	Direct Culture	Not detected
Fungi	Direct Culture	Not detected
Mycoplasma	PCR	Not detected

Cell Type Characterization

Gene	Cell Type	Specification	Results	
NKX6-1	Pancreatic	Positive	Pass	<div style="text-align: center;"> <p style="margin: 0;">NKX6-1 GAPDH</p> <p style="margin: 0;">Size 179bp 208bp</p> </div>
GAPDH	House-Keeping	Positive	Pass	

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.

No. 8, 13520 Crestwood Place
 Richmond BC, Canada V6V 2G2
 T e l : 6 0 4 - 2 4 7 - 2 4 1 6
 F a x : 6 0 4 - 2 4 7 - 2 4 1 4
 w w w . a b m G o o d . c o m