



Certificate of Analysis

Product Description

Product Name	Immortalized Mouse Astrocytes – SV40T (IMA2.1)
Cat Number	T0289
Lot Number	HC7949
Source Organism	Wildtype BALB/c Mice
Source Organ	Cerebral Cortex
Cell Passage Number	21
Viability	>80%
Freezing Medium	Cryopreservation Medium (TM024)
Date of Cryopreservation	April 1, 2014
QC Evaluation Medium	Prigrow III (TM003)

Test for Microbial Contaminants

	Test Method	Results
Bacteria	Direct Culture	Not detected
Fungi	Direct Culture	Not detected
Mycoplasma	PCR	Not detected

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.

1 - 3 6 7 1 Viking Way
 Richmond BC, Canada V4V 2J5
 Tel : 6 0 4 - 2 4 7 - 2 4 1 6
 Fax : 6 0 4 - 2 4 7 - 2 4 1 4
 www.abmGood.com