




Certificate of Analysis

Product Description

Product Name	Luciferase (Luc2) Stable TRAMP Cell Line
Cat Number	T3008
Lot Number	VG0111
Source Organism	Mus musculus
Source Organ	Adenocarcinoma, Prostate
Cell Passage Number	10
Viability	>80%
Freezing Medium	Cryopreservation Medium (TM024)
Date of Cryopreservation	November 12, 2011
QC Evaluation Medium	PriGrow III (TM003)

Test for Microbial Contaminants

	Test Method	Results
Bacteria	Direct Culture	Not detected
Fungi	Direct Culture	Not detected
Mycoplasma	PCR	Not detected



This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.

1-3671 Viking Way
Richmond BC, Canada V4V 2J5
Tel: 604-247-2416
Fax: 604-247-2414
www.abmGood.com