


Certificate of Analysis

Product Description

Product Name	BMP Responsive Immortalized Reporter (BRITER) Cell Line
Cat Number	T3105
Lot Number	AT8054
Source Organism	Bmp2 ^{c/c} ;Bmp4 ^{c/c} ; R26CreER Neonatal Mouse
Source Organ	Calvariae
Cell Passage Number	10
Viability	>80%
Freezing Medium	Cryopreservation Medium (TM024)
Date of Cryopreservation	September 28, 2015
QC Evaluation Medium	Prigrow III (TM003)

Test for Microbial Contaminants

	Test Method	Results
Bacteria	Direct Culture	Not detected
Fungi	Direct Culture	Not detected
Mycoplasma	PCR	Not detected



This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.