


Lenti-PGK-GFP Control Lentivirus



Cat. No.
LV028
Quantity
10ml

Vector	pLenti-III-PGK-GFP Control Vector
Titer	10 ⁶ IU/ml
Selection Marker	Mammalian: Puromycin
Format	Lentivirus
Appearance	Clear Liquid
Storage	1 year when stored at -80°C or lower in a non-frost free freezer.

Caution: This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information (1-866-757-2414).