


## Certificate of Analysis

### Product Description

Product Name	Rat Primary Pulmonary Microvascular Endothelial Cells
Cat Number	T4802
Lot Number	AT8053
Source Organism	<i>Rattus norvegicus</i>
Source Organ	Lung
Cell Passage Number	2
Viability	>80%
Freezing Medium	Cryopreservation Medium (TM024)
Date of Cryopreservation	August 19, 2015
QC Evaluation Medium	Prigrow III (TM003)

### Test for Microbial Contaminants

	Test Method	Results
Bacteria	Direct Culture	Not detected
Fungi	Direct Culture	Not detected
Mycoplasma	PCR	Not detected



This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.