



Certificate of Analysis

Product Description

Product Name	Scrambled AAV siRNA Control Virus (Serotype 2)
Catalog Number	iAAV01502
Lot Number	KL815330
Quantity	5 x 200 µl
Final Storage Buffer	PBS with 5% glycerol

Specifications

	Test Method	Specifications	Results
Viral Titer	qRT-PCR	$> 10^9$ GC/ml	1.17×10^{11} GC/ml

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information.

No. 8, 13520 Crestwood
P l a c e
Richmond BC, Canada V6V
2 G 2
T e l : 6 0 4 2 7 2 1 1 6