

NEW

MegaFi™ Fidelity DNA Polymerase

1000X higher fidelity than regular Taq

Visit: <http://info.abmgood.com/MegaFi>



1000X higher fidelity



Amplifies difficult templates




Effective in long-range PCR




MegaFi™ Fidelity DNA Polymerase has the highest fidelity

A lacZa-containing sequence was amplified using various DNA polymerases, subcloned, and plated on X-gal-containing media. Colonies without PCR amplified lacZa mutations appear dark blue, while those that have PCR introduced frameshift mutations appear white. **abm's** MegaFi™ Fidelity DNA Polymerase is **THE BEST** high fidelity DNA polymerase.




Order the most affordable high fidelity DNA Polymerase:

Product	Size	Cat. No.
MegaFi™ Fidelity DNA Polymerase	400 rxn	G896
MegaFi™ Fidelity 2X PCR MasterMix	800 rxn	G897
MegaFi™ One-Step RT-PCR	100 rxn	G597



Request your FREE sample!



ISO 13485:2016 Certified
Product manufactured under Quality Management System conforming with ISO 13485:2016 as certified by Intertek.