


## Retro-CMV-RFP Retrovirus



**Cat. No.**

RVP005

**Quantity**

1 ml

**Species**

N/A

**Accession Number**

N/A

**Vector**

pRetroG-CMV-RFP-stop Control Vector

**Vector Size**

6960bp

**Insert Size**

N/A

**Sequencing Primer**

5'---CGC AAA TGG GCG GTA GGC GTG---3'

**Tags**

N/A

**Selection Marker**

Mammalian: Puromycin

**Format**

Retrovirus

**Appearance**

Liquid

**Storage Buffer**

Serum-free DMEM

**Storage**

1 year when stored at -80°C.

**Caution:** This product is for research use only and is not intended for therapeutic or diagnostic applications.

Please contact a technical service representative for more information (1-866-757-2414).